

3 steps to success with clariti® 1 day multifocal

Step 1: Start with a new refraction and verification of sensory eye dominance (fogging technique). Convert to spherical equivalent (SE) CL Rx allowing for vertex distance if necessary.

Step 2: Lens selection – modify spherical equivalent (SE) prescription as indicated to select initial trial lens

| | | Spectacle add +0.75 to +1.75 | Spectacle add +2.00 to +2.25 | Spectacle add +2.50 and over |
|------------------------------|------------------|---------------------------------|---------------------------------|---------------------------------|
| | | LOW | | LOW/HIGH |
| Myopes Emmetropes | Dominant Eye | SE | SE | SE +0.25D LOW |
| | Non-dominant Eye | SE | SE +0.50D | SE +0.25D HIGH |
| Hyperopes | Dominant Eye | SE | SE +0.25D | SE +0.25D LOW |
| | Non-dominant Eye | SE +0.25D | SE +0.25D | SE +0.25D HIGH |

SE = Spherical Equivalent LOW = Low add HIGH = High add

Step 3: Allow patients to adapt to lenses for 15 minutes before assessing VA in binocular conditions

- To improve distance VA add +/-0.25D to the dominant eye
- To improve near VA add +0.25D to the non-dominant eye


Example of initial lens selection:

Spectacle prescription, with right eye dominant

RE +3.00D LE +3.00D
Add +1.75


Initial trial lens selection

RE +3.00D LOW LE +3.25D LOW



Clinical tips

- Assess vision in good illumination and with real life scenarios common to wearer e.g. computer, mobile phone, driving distances, etc.
- Do not use phoropter or trial frame when assessing/improving vision. Use handheld trial lenses without occlusion.
- Subjective quality of vision is usually more important than Snellen acuity



| clariti® 1 day multifocal product specifications | |
|--|--|
| Material | somofilcon A |
| Water content | 56% |
| Base curve | 8.6 mm |
| Diameter | 14.1 mm |
| Sphere power range | +5.00D to -6.00D (0.25D steps) |
| Power Additions | LOW: spec add up to +2.25D HIGH: spec add +2.50D and over |
| Center thickness (@ -3.00D) | 0.07 mm |
| Dk | 60 |
| Dk/t (@ -3.00D) | 86 |
| Modulus | 0.5 MPa |
| UV filter* inhibitor | Yes |
| Pack size | 30-pack / 90-pack 10-pack trial |
| Lens design | Center-near sphere |

*Warning: UV-absorbing contact lenses are not substitutes for protective UV-absorbing eyewear, such as UV-absorbing goggles or sunglasses, because they do not completely cover the eye and surrounding area. Patients should continue to use UV-absorbing eyewear as directed.