

Product Reference Guide

Look to the people, products, and services of CooperVision®.

CONTENTS

Spheres	2 – 3
Torics	4 – 6
Multifocals	7 – 8
Toric Multifocals	8
Return & Exchange Policies	9

Customer Service: 800 341 2020
www.coopervision.com



CooperVision®
Live Brightly.®

Product Reference Guide



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS



OPTIMIZED BALLAST DESIGN













UV BLOCKING



DESIGNED FOR MYOPIA CONTROL

SPHERE LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
1-Day Disposables										
clariti® 1 day 	somofilcon A / 56%	8.6	14.1	-0.50D to -10.00D +0.50D to +8.00D (0.50D steps after +/- 6.00D) No Plano	SiHy	Daily Disposable	86	90-pack and 30-pack	56% 1 day Silicone Hydrogel sphere / 10 pk	
MyDay® 	stenfilcon A / 54%	8.4	14.2	+8.00D to -12.00D (0.50D steps after +5.00D and -6.00D) No Plano	SiHy	Daily Disposable	100	90-pack and 180-pack	54% 1 day asphere / 5pk	
MyDay Energys® 	stenfilcon A / 54%	8.4	14.2	+8.00D to -12.00D (0.50D steps after +5.00D and -6.00D) No Plano	SiHy	Daily Disposable	100	90-pack	54% 1 day E asphere / 10pk	
Proclear® 1 day 	omafilcon A / 60%	8.7	14.2	+8.00D to -12.00D (0.50D steps after +5.00D and -6.00D) No Plano	2	Daily Disposable	25	90-pack	60% 1 day sphere / 10pk	
MiSight® 1 day 	omafilcon A / 60%	8.7	14.2	-0.50D to -7.00D (0.50D steps after -6.00D) No plano	2	Daily Disposable	25	90-pack	Myopia Management ActivControl™ Technology / 10pk	

Product Reference Guide



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS



OPTIMIZED BALLAST DESIGN



UV BLOCKING



DESIGNED FOR MYOPIA CONTROL

SPHERE LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
2-Week Replacement Lenses										
 Avaira Vitality™	fanfilcon A / 55%	8.4	14.2	+8.00D to -12.00D (0.50D steps after +/-6.00D) No plano	SiHy	Daily wear	110	6-pack	fanfilcon A / single	
 Biomedics® 55 Premier asphere	ocufilcon D / 55%	8.6	14.2	-0.25D to -10.00D (0.50D steps after -6.00D) No plano	4	Daily wear	28.1	6-pack	ocufilcon D 55% asphere / 3pk	
		8.8	14.2	+0.25D to +8.00D (0.50D steps after +5.00D) No plano						
Monthly Replacement Lenses										
 Biofinity®	comfilcon A / 48%	8.6	14.0	+8.00D to -12.00D (0.50D steps after ±6.00D) No plano	SiHy	Daily Wear or Extended Wear up to 6 nights/ 7 days; monthly replacement	160	6-pack	comfilcon A / single	
 Biofinity® XR	comfilcon A / 48%	8.6	14.0	+8.50D to +15.00D -12.50D to -20.00D (0.50D steps) No plano	SiHy	Daily Wear or Extended Wear up to 6 nights/ 7 days; monthly replacement	160	6-pack	comfilcon A XR / single	
 Biofinity Energys®	comfilcon A / 48%	8.6	14.0	+8.00D to -12.00D (0.50D steps after ±6.00D) No plano	SiHy	Daily Wear or Extended Wear up to 6 nights/ 7 days; monthly replacement	160	6-pack	Biofinity Energys® / single	
 Proclear® sphere	omafilcon B / 62%	8.6	14.2	+20.00D to -20.00D (0.50D steps after +/-6.00D)	2	Daily wear	36	6-pack	Proclear / 3pk	
		8.2	14.2	Plano to -10.00D (0.50D steps after -6.00D)						

Product Reference Guide



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS



OPTIMIZED BALLAST DESIGN







UV BLOCKING



DESIGNED FOR MYOPIA CONTROL

TORIC LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
1-Day Disposables												
clariti® 1 day toric 	somofilcon A / 56%	8.6	14.3	plano to -6.00D (0.25D steps)	-0.75, -1.25, -1.75	Full circle in 10 degrees	SiHy	Daily Disposable	57	90-pack and 30-pack	56% 1 day Silicone Hydrogel toric / 5pk	
				-6.50D to -9.00D (0.50D steps)	-2.25	10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°						
				+0.25D to +4.00D (0.25D steps)	-0.75, -1.25, -1.75	10°, 20°, 60°, 70°, 80°, 90°, 100°, 110°, 120°, 160°, 170°, 180°						
					-2.25	10°, 20°, 90°, 160°, 170°, 180°						
MyDay® toric 	stenfilcon A / 54%	8.6	14.5	+8.00D to -10.00D (0.50D steps after +/-6.00D)	-0.75, -1.25, -1.75, -2.25	Full circle in 10° steps	SiHy	Daily Disposable	80	90-pack	54% 1 day/ Silicone Hydrogel toric / 5pk	

Product Reference Guide



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS



OPTIMIZED BALLAST DESIGN














UV BLOCKING



DESIGNED FOR MYOPIA CONTROL

TORIC LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
2-Week Replacement Lenses												
Avaira Vitality™ toric 	fanfilcon A / 55%	8.5	14.5	+8.00D to -10.00D (0.50D steps after ± 6.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	SiHy	Daily wear	90	6-pack	fanfilcon A / single	     
Biomedics® Toric 	ocufilcon D / 55%	8.7	14.5	+6.00D to -10.00D (0.50D steps after -6.00D and +5.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	4	Daily wear	17.8	6-pack	ocufilcon D toric / single	  

Product Reference Guide



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS



OPTIMIZED BALLAST DESIGN





















UV BLOCKING



DESIGNED FOR MYOPIA CONTROL

TORIC LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/ Pack Size	Features
Monthly Replacement Lenses												
Biofinity® toric 	comfilcon A / 48%	8.7	14.5	+8.00D to -10.00D (0.50D steps after +/-6.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	SiHy	Daily Wear or Extended Wear up to 6 nights/ 7 days; monthly replacement	116	6-pack	comfilcon A toric / single	    
Biofinity® XR toric 	comfilcon A / 48%	8.7	14.5	+8.50D to +20.00D and -10.50D to -20.00D (0.50D steps) +20.00D to -20.00D (0.50D steps after +/-6.00D)	-0.75 -1.25 -1.75 -2.25 -2.75 -3.25 -3.75 -4.25 -4.75 -5.25 -5.75	5° to 180° (in 5° steps)	SiHy	Daily Wear or Extended Wear up to 6 nights/ 7 days; monthly replacement	116	6-pack	comfilcon A toric / single	    
Proclear® toric 	omafilcon B / 62%	8.8 8.4	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)*	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	2	Daily wear	25	6-pack	Proclear® toric / single	 
Proclear® toric XR 	omafilcon B / 62%	8.8 8.4	14.4	+6.25D to +10.00D (0.50D steps after +6.50D) -8.50D to -10.00D (0.50D steps)* +10.00D to -10.00D (0.50D steps after +/-6.50D)*	-0.75 -1.25 -1.75 -2.25 -2.75 to -5.75 (0.50D steps)	5° to 180° (5° steps)	2	Daily wear	39	6-pack	Proclear® toric / single	 

Product Reference Guide



VISIBILITY
TINT



LOW
MODULUS



NATURALLY
WETTABLE



HIGH OXYGEN
TRANSMISSIBILITY



ASPHERIC
OPTICS



OPTIMIZED
BALLAST DESIGN


























UV
BLOCKING



DESIGNED FOR
MYOPIA CONTROL

MULTIFOCAL LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	ADD Power	Design	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
1-Day Disposables												
clariti® 1 day multifocal 	somofilcon A / 56%	8.6	14.1	+5.00D to -6.00D in 0.25 D steps	Low: up to +2.25D High: +2.50D to +3.00D	center near, progressive intermediate and peripheral distance	SiHy	Daily Disposable	86	90-pack and 30-pack	56% 1 day Silicone Hydrogel multifocal / 10 pk	    
MyDay® multifocal 	stenfilcon A / 54%	8.4	14.2	+8.00D to -10.00D (0.25 D steps) -10.50D to -12.00D (0.50 D steps)	Low: +0.75D to +1.25D Mid: +1.50D to +1.75D High: +2.00D to +2.50D	Binocular Progressive System™	SiHy	Daily Disposable	100	90-pack	54% 1 day Silicone Hydrogel multifocal / 10pk	    
Proclear® 1 day multifocal 	omafilcon A / 60%	8.7	14.2	+6.00D to -10.00D (0.50D steps after -6.00D)	SINGLE POWER PROFILE	center near aspheric	2	Daily Disposable	28	90-pack	60% 1 day PC multifocal/ 10 pk	 
Monthly Replacement Lenses												
Biofinity® multifocal 	comfilcon A / 48%	8.6	14.0	+6.00D to -10.00D (0.50D steps after -6.00D)	+1.00 +1.50 +2.00 +2.50	D lens N lens	SiHy	Daily Wear or Extended Wear up to 6 nights/ 7 days; monthly replacement	142	6-pack	comfilcon A multifocal / single	   
Proclear® multifocal 	omafilcon B / 62%	8.7	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)	+1.00 +1.50 +2.00 +2.50	D lens N lens	2	Daily wear	17	6-pack	Proclear® multifocal / single	 

Product Reference Guide



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS



OPTIMIZED BALLAST DESIGN






UV BLOCKING












DESIGNED FOR MYOPIA CONTROL

MULTIFOCAL LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	ADD Power	Design	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Monthly Replacement Lenses												
Proclear® multifocal XR 	omafilcon B / 62%	8.7	14.4	+6.25D to +20.00D -8.50D to -20.00D (0.50D steps after +/- 6.50D)	+1.00 +1.50 +2.00 +2.50	D lens N lens	2	Daily wear	17	6-pack	Proclear multifocal / single	 
		8.7	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D)	+3.00 +3.50 +4.00							
		8.4	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D)	+1.00, +1.50, +2.00, +2.50, +3.00, +3.50, +4.00							

TORIC MULTIFOCAL LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	ADD Power	Design	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Monthly Replacement Lenses														
Biofinity® toric multifocal 	comfilcon A / 48%	8.7	14.5	+10.00D to -10.00D (0.50D steps after +/- 6.00D)	-0.75D to -5.75D (0.50D steps)	5°-180°, 5° steps	+1.00 +1.50 +2.00 +2.50	D lens N lens	SiHy	Daily Wear or Extended Wear up to 6 nights/7 days; monthly replacement	116	6-pack	Biofinity® toric multifocal / single	    
Proclear® multifocal toric 	omafilcon B / 62%	8.8 8.4	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D)	-0.75 to -5.75 (0.50D steps)	5° to 180° (5° steps)	+1.00 to +4.00 (0.50D steps)	D lens N lens	2	Daily wear	38.6	6-pack	Proclear® multifocal toric / single	 

Return and Exchange Policies

SIMPLY SATISFIED GUARANTEE™

CooperVision's Simply Satisfied Guarantee™ allows the return of unopened boxes of CooperVision lenses purchased directly by customers from CooperVision.



Brands eligible for the Simply Satisfied Guarantee™

Avaira Vitality® **Biofinity®** **MyDay®** **MiSight® 1 day**
clariti® 1 day **Biomedics®** **Proclear®†**

† *Proclear® toric XR and Proclear® multifocal XR are not eligible for return or exchange.*

All returns must meet the requirements below:

- Product must be returned within 100 days of the original invoice date.
 - Discontinued product must be returned 30 days prior to product discontinuance date.
- A copy of the entire original corresponding CooperVision invoice must be included and indicate the reason for return.
- Product must have a minimum of 1 year dating of the expiration date shown on the packaging.
- Product must be unopened, intact and in resalable condition (as determined by CooperVision), free of writing, markings of any kind, and stickers not applied by CooperVision.
 - Lenses in packages with the end flaps removed, due to patient rebate redemption or any other reason are not eligible for return.
- Lenses must be received in the original packaging, including the original blister or vial.

* *CooperVision may destroy product that does not meet the return requirements without credit or refund.*

MISIGHT® 1 DAY OPENED BOX EXCHANGE POLICY

Allows for exchange of opened boxes of MiSight® 1 day within 100 days of the original invoice date. Product must have a minimum of 1 year dating of the expiration date shown on the packaging and there must be at least 75 unopened, intact lens blisters in each 90 pack to be eligible. A copy of the original corresponding CooperVision invoice and reason for exchange are required. Packages with the end flaps removed, due to patient rebate redemption or any other reason are not eligible for exchange.

* *CooperVision may destroy product that does not meet the exchange requirements without credit or refund.*

UNOPENED EXCHANGE POLICY

Our Unopened Exchange Policy allows exchanges of most lenses currently offered for sale and purchased directly from CooperVision.

Products eligible for Simply Satisfied Guarantee™ are eligible for an unopened exchange.

† *Proclear® toric XR and Proclear® multifocal XR are not eligible for return or exchange.*

All exchanges must meet the requirements below:

- Product must be exchanged for the same lens type (i.e. Avaira Vitality™ 6-pack for Avaira Vitality™ 6-pack, Biofinity Energys® 6-pack for Biofinity Energys® 6-pack.)
- The same lens type must have been ordered within the last 100 days of the original invoice date.
- A copy of the entire original corresponding CooperVision invoice must be included and indicate the reason for exchange.
- Product must have a minimum of 1 year dating of the expiration date shown on the packaging.
- Product must be unopened, intact and in resalable condition (as determined by CooperVision), free of writing, markings of any kind, and stickers not applied by CooperVision.
 - Lenses in packages with the end flaps removed, due to patient rebate redemption or any other reason are not eligible for return.

* *CooperVision may destroy product that does not meet the exchange requirements without credit or refund.*

MADE-TO-ORDER LENSES:

- CooperVision offers some lenses that are made-to-order and will ship within 10 days from the time of order. These lenses are Proclear® toric XR, Proclear® multifocal XR and Proclear® multifocal toric. Biofinity® XR toric* will ship within 14-21 days from the time of order. (*Biofinity® XR toric is an MTO product with limited Rx in stock).

UV NOTICE:

- UV-absorbing contact lenses are NOT substitutes for protective UV-absorbing eyewear, such as UV-absorbing goggles or sunglasses, as they do not completely cover the eye and surrounding area. Patients should continue to use UV-absorbing eyewear as directed.