

Product Reference Guide

Look to the people, products, and services of CooperVision®.

CONTENTS

Spheres	2 – 5
Torics	6 – 8
Multifocals	9 – 10

Customer Service: 800 341 2020
www.coopervision.com



CooperVision®
Live Brightly.®



VISIBILITY
TINT



LOW
MODULUS



NATURALLY
WETTABLE



HIGH OXYGEN
TRANSMISSIBILITY



ASPHERIC
OPTICS
























OPTIMIZED
BALLAST DESIGN



UV
BLOCKING

SPHERE LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
1-Day Disposables										
clariti® 1 day 	somofilcon A / 56%	8.6	14.1	-0.50D to -10.00D +0.50D to +8.00D (0.50D steps after ±6.00D) No Plano	SiHy	Daily Disposable	86	90-pack	1 day silicone hydrogel/10 pk	    
ClearSight™ 1-Day 	ocufilcon B/ 52%	8.7	14.2	+6.00D to -10.00D (0.50D steps after +5.00D and -6.00D) No Plano	4	Daily Disposable	24	90-pack and 30-pack	52% 1-Day Trial lenses (ocufilcon B)/ 10pk	  
MyDay® 	stenfilcon A/ 54%	8.4	14.2	+8.00D to -12.00D (0.50D steps after +5.00D and -6.00D) No Plano	SiHy	Daily Disposable	100	90-pack and 180-pack	54% 1 Day Silicone Hydrogel/10pk	     
Proclear® 1 day 	omafilcon A / 60%	8.7	14.2	+8.00D to -12.00D (0.50D steps after +5.00D and -6.00D) No Plano	2	Daily Disposable	28	90-pack and 30-pack	60% 1-Day Trial lenses (omafilcon A)/ 10pk	  



VISIBILITY
TINT



LOW
MODULUS



NATURALLY
WETTABLE



HIGH OXYGEN
TRANSMISSIBILITY



ASPHERIC
OPTICS






















OPTIMIZED
BALLAST DESIGN



UV
BLOCKING

SPHERE LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
2-Week Replacement Lenses										
Avaira®  Product being discontinued 5/1/18	enfilcon A / 46%	8.4	14.2	+0.25D to +8.00D (0.50D steps after +6.00D)	SiHy	Daily	125	6-pack	No Trials	     
		8.5	14.2	-0.25D to -12.00D (0.50D steps after -6.00D)						
Avaira Vitality™ 	fanfilcon A / 55%	8.4	14.2	+8.00D to -12.00D (0.50D steps after +/-6.00D) No plano	SiHy	Daily	110	6-pack	fanfilcon A/ single	     
Biomedics® 55 Premier asphere 	ocufilcon D/ 55%	8.6	14.2	-0.25D to -10.00D (0.50D steps after -6.00D)	4	Daily or Extended	28.1	6-pack	55 UV aspheric contact lens/ 6	   
		8.9		8.8						



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS














OPTIMIZED BALLAST DESIGN



UV BLOCKING

SPHERE LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Monthly Replacement Lenses										
Biofinity® 	comfilcon A/ 48%	8.6	14.0	+8.00D to -12.00D (0.50D steps after +/-6.00D)	SiHy	Daily or 6 nights/ 7 days Extended	160	6-pack	Biofinity (comfilcon A)/ single	
Biofinity® XR 	comfilcon A/ 48%	8.6	14.0	+8.50D to +15.00D -12.50D to -20.00D (0.50D steps)	SiHy	Daily or 6 nights/ 7 days Extended	160	6-pack	Biofinity (comfilcon A)/ single	
Biofinity® Energys™ 	comfilcon A / 48%	8.6	14.0	+8.00D to -12.00D (0.50D steps after +/-6.00D)	SiHy	Daily or 6 nights/ 7 days Extended	160	6-pack	comfilcon A asphere/ single	
Expressions® Colors 	methafilcon A/ 55%	8.7	14.4	+4.00D to -6.00D (0.25D steps)	4	Daily	37.6	6-pack	No Trials	Colors available: Jade, Blue Topaz, Aqua, Blue, Green, Gray, Hazel, Brown
Frequency® 55 Aspheric 	methafilcon A/ 55%	8.7	14.4	+8.00D to -10.00D (0.50D steps after -8.00D)	4	Daily	20.9	6-pack	Aspheric (methafilcon A)/ 3	
		8.4	14.4	-0.25D to -8.00D (0.25D steps)						
Proclear® sphere 	omafilcon B / 62%	8.6	14.2	+20.00D to -20.00D (0.50D steps after +/-6.00D)	2	Daily	52.3	6-pack	Proclear (omafilcon A)/ 3	
		8.2	14.2	Plano to -10.00D (0.50D steps after -6.00D)						



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS



OPTIMIZED BALLAST DESIGN



UV BLOCKING

SPHERE LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Quarterly or Conventional Replacement Lenses										
Preference® <i>Product being discontinued 11/1/18</i>	tetrafilcon A/ 43%	8.4 8.7	14.4	-0.25D to -10.00D (0.50D steps after -6.50D)	1	Daily	26.6	4-pack vials	No Trials	
Preference® Standard <i>Product being discontinued 11/1/18</i>	tetrafilcon A/ 43%	8.3 8.6	14.0 14.0	+6.00D to -10.00D (0.50D steps after -6.50D)	1	Daily	20.9	4-pack vials	No Trials	
Hydrasoft® Sphere	methafilcon B/ 55%	8.3 8.6	14.2	+10.00D to -30.00D (0.25D steps)	4	Daily	18.8	Single and 3-pack vials	No Trials (Singles available for purchase)	
		8.3 8.6 8.9 9.2	15.0				15.7			
Hydrasoft® Sphere Thin	methafilcon B/ 55%	8.3 8.6	14.2	+10.00D to -30.00D (0.25D steps)	4	Daily	31.3	Single and 3-pack vials	No Trials (Singles available for purchase)	
		8.3 8.6 8.9 9.2	15.0							
Hydrasoft® Aphakic	methafilcon B/ 55%	8.3 8.6	14.2	+10.25D to +30.00D (0.25D steps)	4	Daily		Single and 3-pack vials	No Trials (Singles available for purchase)	
		8.3 8.6 8.9 9.2	15.0							
Hydrasoft® Aphakic Thin	methafilcon B/ 55%	8.3 8.6	14.2	+10.25D to +30.00D (0.25D steps)	4	Daily		Single and 3-pack vials	No Trials (Singles available for purchase)	
		8.3 8.6 8.9 9.2	15.0							



VISIBILITY
TINT



LOW
MODULUS



NATURALLY
WETTABLE



HIGH OXYGEN
TRANSMISSIBILITY



ASPHERIC
OPTICS




































OPTIMIZED
BALLAST DESIGN



UV
BLOCKING

TORIC LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
1-Day Disposables												
clariti® 1 day toric 	somofilcon A / 56%	8.6	14.3	Plano to - 9.00D (0.50D steps after -6.00D)	-0.75, -1.25, -1.75	10°, 20°, 60° 70°, 80°, 90° 100°, 110°, 120° 160°, 170°, 180°	SiHy	Daily Disposable	57	90-pack and 30-pack	1 day toric silicone hydrogel I/ 5pk	     
				+0.25D to +4.00D (0.25D steps)		20°, 70°, 90° 110°, 160°, 180°						
				Plano to -9.00D (0.50D steps after -6.00D)	-2.25	10°, 20°, 90° 160°, 170°, 180°						
MyDay® toric 	stenfilcon A/ 54%	8.6	14.5	Plano to -10.00D (0.50D steps after -6.00D)	-0.75 -1.25 -1.75	160°, 170°, 180°, 10°, 20°, 70°, 80°, 90°, 100°, 110°	SiHy	Daily Disposable	80	90-pack	54% 1 day/ Silicone Hydrogel toric/ 5pk	      
				Coming Soon: parameters extending to +6.00D		-2.25						
2-Week Replacement Lenses												
Avaira® toric  Product being discontinued 5/1/18	enfilcon A / 46%	8.5	14.5	+6.00D to -10.00D (0.50 steps after +/-6.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	SiHy	Daily	91	6-pack	No Trials	     
Avaira Vitality™ toric 	fanfilcon A/ 55%	8.5	14.5	+8.00D to -10.00D (0.50D steps after ± 6.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	SiHy	Daily	90	6-pack	fanfilcon A/ single	     
Biomedics® Toric 	ocufilcon D/ 55%	8.7	14.5	+6.00D to -9.00D (0.50 steps after -6.00D and +5.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	4	Daily or Extended	17.8	6-pack	toric contact lens/ single	  



VISIBILITY
TINT



LOW
MODULUS



NATURALLY
WETTABLE



HIGH OXYGEN
TRANSMISSIBILITY



ASPHERIC
OPTICS













OPTIMIZED
BALLAST DESIGN



UV
BLOCKING

TORIC LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Vertex Toric® 	methafilcon A/ 55%	8.6	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	4	Daily	17.9	6-pack	Vertex toric/ single	
Monthly Replacement Lenses												
Biofinity® toric 	comfilcon A/ 48%	8.7	14.5	+8.00D to -10.00D (0.50D steps after +/-6.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	SiHy	Daily or up to 7 days extended	116	6-pack	comfilcon A toric / single	
Biofinity® XR toric 	comfilcon A/ 48%	8.7	14.5	+10.00D to +8.50D (0.50D steps) +10.00D to -10.00D (0.50D steps after +/- 6.00D)	-0.75 -1.25 -1.75 -2.25 -2.75 -3.25 -3.75 -4.25 -4.75 -5.25 -5.75	5° to 180° (in 5° steps)	SiHy	Daily or up to 7 days extended	116	6-pack	comfilcon A toric / single	
Frequency® 55 Toric  <i>Product being discontinued 11/1/2018</i>	methafilcon B/ 55%	8.4 8.7	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	4	Daily	15.7	6-pack	No Trials	
Frequency® 55 Toric XR  <i>Product being discontinued 11/1/2018</i>	methafilcon B/ 55%	8.4 8.7	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)	-2.75 to -5.75 (0.50D steps)	5° to 180° (5° steps)	4	Daily	15.7	6-pack	No Trials	



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS















OPTIMIZED BALLAST DESIGN



UV BLOCKING

TORIC LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Proclear® toric 	omafilcon B/ 62%	8.8 8.4	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)*	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	2	Daily	30.9	6-pack	Proclear toric/ single	 
Proclear® toric XR 	omafilcon B/ 62%	8.8 8.4	14.4	+6.25D to +10.00D (0.50D steps after +6.50D) -8.50D to -10.00D (0.50D steps)*	-0.75 -1.25 -1.75 -2.25	5° to 180° (5° steps)	2	Daily	30.9	6-pack	Proclear toric/ single	 
		8.8 8.4	14.4	+10.00D to -10.00D (0.50D steps after +/-6.50D)*	-2.75 to -5.75 (0.50D steps)							
Quarterly or Conventional Replacement Lenses												
Hydrasoft® Toric	methafilcon B/ 55%	8.3 8.6	14.2	+30.00D to -30.00D (0.25D steps)	-0.50D to -10.00D (0.25D steps)	1° to 180° (1° steps)	4	Daily	14.5	Single and 3-pack vials	No Trials (Singles available for purchase)	
		8.3 8.6 8.9 9.2	15.0									
Hydrasoft® Toric Thin	methafilcon B/ 55%	8.3 8.6	14.2	+30.00D to -30.00D (0.25D steps)	-0.50D to -10.00D (0.25D steps)	1° to 180° (1° steps)	4	Daily	14.5	Single and 3-pack vials	No Trials (Singles available for purchase)	
		8.3 8.6 8.9 9.2	15.0									
Preference® Toric <i>Product being discontinued 11/1/2018</i>	tetrafilcon A/ 43%	8.4 8.7	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)	-0.75 -1.25 -1.75 -2.25	5° to 180° (5° steps)	1	Daily	10.3	4-pack vials	No Trials	 
Preference® Toric XR <i>Product being discontinued 11/1/2018</i>	tetrafilcon A/ 43%	8.4 8.7	14.4	+6.00D to -9.50D (0.50 steps after -6.50D)**	-2.75 to -9.75 (0.25 steps)	5° to 180° (5° steps)	1	Daily	10.3	4-pack vials	Preference toric / single <i>Trials being discontinued 5/1/2018</i>	 

* Proclear toric XR is available in the full range of sphere powers and cylinder powers. Sphere powers over +6.00D and -8.00D would be considered XR for all listed cylinder powers. Lenses with a -2.75 to -5.75 cylinder with sphere powers UNDER +6.00D or -8.00D also would be within the XR range.

** Preference Toric XR is available in the full range of sphere powers and cylinder powers. Sphere powers over -8.00D would be considered XR for all listed cylinder powers. Lenses with a -2.75 to -9.75 cylinder with sphere powers UNDER +6.00D or -8.00D also would be within the XR range.



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS




















OPTIMIZED BALLAST DESIGN



UV BLOCKING

MULTIFOCAL LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	ADD Power	Design	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
1-Day Disposables												
clariti® 1 day multifocal 	somofilcon A / 56%	8.6	14.1	+5.00D to -6.00D (0.25D steps)	Low: up to +2.25D High: +2.50D to +3.00D	center near, progressive intermediate and peripheral distance	SiHy	Daily Disposable	86	90-pack and 30-pack	1 day multifocal silicone hydrogel / 10 pk	    
Proclear® 1 day multifocal 	omafilcon A / 60%	8.7	14.2	+6.00D to -10.00D (0.50D steps after -6.00D)	SINGLE POWER PROFILE	center near aspheric	2	Daily Disposable	28	90-pack and 30-pack	60% 1 Day Proclear multifocal / 10 pk	 
Monthly Replacement Lenses												
Biofinity® multifocal 	comfilcon A / 48%	8.6	14.0	+6.00D to -10.00D (0.50D steps after -6.00D)	+1.00, +1.50, +2.00, +2.50	D lens N lens	SiHy	Daily or up to 7 days extended	142	6-pack	comfilcon A multifocal / single	   
Proclear® multifocal 	omafilcon B / 62%	8.7	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)	+1.00, +1.50, +2.00, +2.50	D lens N lens	2	Daily	21.3	6-pack	Proclear multifocal / single	 

Product Reference Guide



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS







OPTIMIZED BALLAST DESIGN



UV BLOCKING

MULTIFOCAL LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	ADD Power	Design	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Proclear® multifocal XR 	omafilcon B/ 62%	8.7	14.4	+6.25D to +20.00D -8.50D to -20.00D (0.50D steps after +/- 6.50D)	+1.00, +1.50, +2.00, +2.50	D lens N lens	2	Daily	21.3	6-pack	Proclear multifocal / single	
		8.7	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D)	+3.00, +3.50, +4.00							
		8.4	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D)	+1.00, +1.50, +2.00, +2.50, +3.00, +3.50, +4.00							
Proclear® multifocal toric 	omafilcon B/ 62%	8.8 8.4	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D) CYL: -0.75 to -5.75 (0.50D steps) AXIS: 5° to 180° (5° steps)	+1.00 to +4.00 (0.50 steps)	D lens N lens	2	Daily	48.6	6-pack	Proclear multifocal toric / single	

"IT'S OKAY"SM GUARANTEE:



CooperVision's "It's Okay"SM Guarantee allows the return of CooperVision lenses purchased directly from CooperVision.

Brands eligible for the "It's Okay"SM Guarantee"

Avaira Vitality™ Brand	Biofinity® Brand	Biomedics® Brand
clariti® Brand	ClearSight™ Brand	MyDay® Brand
Proclear® Brand	Hydrasoft® Brand	Vertex® Brand
	<small>(Single vials only)</small>	

[†] Hydrasoft 3-pack is not eligible for return or exchange.

For brands not listed above, please refer to the Unopened Exchange Policy.

All returns must meet the requirements below:

- Product must be returned within 100 days of the original invoice date.
 - Discontinued product must be returned within 30 days of product discontinuance date.
- A copy of the original corresponding CooperVision invoice must be included and indicate the reason for return.
- Product must have a minimum of 1 year dating of the expiration date shown on the packaging.
- Unopened product must be intact and in resalable condition, free of writing, stickers or markings of any kind.
 - Lenses in packages with the end flaps removed, due to patient rebate redemption or any other reason are not eligible for return.
- Multi packs must contain a minimum of lenses as indicated below.

6 Pack: 5 • 30 Pack: 25 • 90 Pack: 85 • 180 Pack: 170

** CooperVision may destroy product not eligible for return.*

Pack Size	Minimum Number of Unopened Lenses Required
6 Pack	5
30 Pack	25
90 Pack	85

- All lenses must be returned with an invoice or completed return form matching the product returned and indicate a reason for return.

UNOPENED EXCHANGE POLICY

Our Unopened Exchange Policy allows exchanges of most lenses currently offered for sale and purchased directly from CooperVision.

Products eligible for "It's Okay"SM Guarantee and the additional products listed below are eligible for an unopened exchange.

Frequency® 55 Brand Expressions® Brand Preference® Brand

[†] Hydrasoft 3-pack is not eligible for exchange or return.

All exchanges must meet the requirements below:

- Product must be exchanged for the same lens type (i.e. Avaira Vitality™ 6 pack for Avaira Vitality™ 6 pack, Biofinity Energys™ 6 pack for Biofinity Energys™ 6 pack.)
- The same lens type must have been ordered within the last 100 days.
- Product must be unopened, intact and in resalable condition, free of writing, stickers or markings of any kind.
 - Lenses in packages with the end flaps removed, due to patient rebate redemption or any other reason are not eligible for return.
- Product must have a minimum of 1 year dating of the expiration date shown on the packaging.

** CooperVision may destroy product not eligible for exchange.*

MADE-TO-ORDER LENSES:

CooperVision offers some lenses that are made-to-order and will ship within 10 days from the time of order. These lenses are Proclear® toric XR, Proclear® multifocal XR, Proclear® multifocal toric, Frequency® 55 Toric XR, and Preference, Preference® toric, Preference® toric XR, and Biofinity® XR toric. Hydrasoft® and Hydrasoft® Toric will ship within 14-21 days from the time of order.

UV NOTICE:

UV-absorbing contact lenses are NOT substitutes for protective UV-absorbing eyewear, such as UV-absorbing goggles or sunglasses, as they do not completely cover the eye and surrounding area. Patients should continue to use UV-absorbing eyewear as directed.

