

# Product Reference Guide

Look to the people, products, and services of CooperVision®.

## CONTENTS

Spheres .....	2 – 5
Torics .....	6 – 8
Multifocals .....	9 – 10

Customer Service: 800 341 2020  
[www.coopervision.com](http://www.coopervision.com)



CooperVision®  
Live Brightly.®



VISIBILITY  
TINT



LOW  
MODULUS



NATURALLY  
WETTABLE



HIGH OXYGEN  
TRANSMISSIBILITY



ASPHERIC  
OPTICS
























OPTIMIZED  
BALLAST DESIGN



UV  
BLOCKING

SPHERE LENSES

Brand	Material/ H <sub>2</sub> O Content	Base Curve (mm)	Dia (mm)	Sphere Power	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
<b>1-Day Disposables</b>										
clariti® 1 day 	somofilcon A / 56%	8.6	14.1	-0.50D to -10.00D +0.50D to +8.00D (0.50D steps after ±6.00D) No Plano	SiHy	Daily Disposable	86	90-pack and 30-pack	1 day silicone hydrogel/10 pk	    
ClearSight™ 1-Day 	ocufilcon B / 52%	8.7	14.2	+6.00D to -10.00D (0.50D steps after +5.00D and -6.00D) No Plano	4	Daily Disposable	24	90-pack	52% 1-Day Trial lenses (ocufilcon B)/ 10pk	  
MyDay® 	stenfilcon A / 54%	8.4	14.2	+8.00D to -12.00D (0.50D steps after +5.00D and -6.00D) No Plano	SiHy	Daily Disposable	100	90-pack and 180-pack	54% 1 Day Silicone Hydrogel/10pk	     
Proclear® 1 day 	omafilcon A / 60%	8.7	14.2	+8.00D to -12.00D (0.50D steps after +5.00D and -6.00D) No Plano	2	Daily Disposable	28	90-pack	60% 1-Day Trial lenses (omafilcon A)/ 10pk	  



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS















OPTIMIZED BALLAST DESIGN



UV BLOCKING

SPHERE LENSES

Brand	Material/ H <sub>2</sub> O Content	Base Curve (mm)	Dia (mm)	Sphere Power	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
2-Week Replacement Lenses										
Avaira Vitality™ 	fanfilcon A / 55%	8.4	14.2	+8.00D to -12.00D (0.50D steps after +/-6.00D) No plano	SiHy	Daily	110	6-pack	fanfilcon A/ single	     
Biomedics® 55 Premier asphere 	ocufilcon D / 55%	8.6 8.9	14.2	-0.25D to -10.00D (0.50D steps after -6.00D)	4	Daily or Extended	28.1	6-pack	55 UV aspheric contact lens/ 3	   
		8.8	14.2	+0.25D to +6.00D (0.50D steps after +5.00D)						



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS














OPTIMIZED BALLAST DESIGN



UV BLOCKING

SPHERE LENSES

Brand	Material/ H <sub>2</sub> O Content	Base Curve (mm)	Dia (mm)	Sphere Power	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Monthly Replacement Lenses										
 Biofinity®	comfilcon A / 48%	8.6	14.0	+8.00D to -12.00D (0.50D steps after +/-6.00D)	SiHy	Daily or 6 nights/ 7 days Extended	160	6-pack	comfilcon A/ single	
 Biofinity® XR	comfilcon A / 48%	8.6	14.0	+8.50D to +15.00D -12.50D to -20.00D (0.50D steps)	SiHy	Daily or 6 nights/ 7 days Extended	160	6-pack	comfilcon A/ single	
 Biofinity Energys®	comfilcon A / 48%	8.6	14.0	+8.00D to -12.00D (0.50D steps after +/-6.00D)	SiHy	Daily or 6 nights/ 7 days Extended	160	6-pack	Biofinity Energys® / single	
 Expressions® Colors <i>Product being discontinued 11/1/19</i>	methafilcon A / 55%	8.7	14.4	+4.00D to -6.00D (0.25D steps)	4	Daily	37.6	6-pack	No Trials	Colors available: Jade, Blue Topaz, Aqua, Blue, Green, Gray, Hazel, Brown
 Frequency® 55 Aspheric <i>Product being discontinued 11/1/19</i>	methafilcon A / 55%	8.7	14.4	+8.00D to -10.00D (0.50D steps after -8.00D)	4	Daily	20.9	6-pack	No Trials	
		8.4	14.4	-0.25D to -8.00D (0.25D steps)						
 Proclear® sphere	omafilcon B / 62%	8.6	14.2	+20.00D to -20.00D (0.50D steps after +/-6.00D)	2	Daily	52.3	6-pack	Proclear (omafilcon A)/ 3	
		8.2	14.2	Plano to -10.00D (0.50D steps after -6.00D)						



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS



OPTIMIZED BALLAST DESIGN



UV BLOCKING

SPHERE LENSES

Brand	Material/ H <sub>2</sub> O Content	Base Curve (mm)	Dia (mm)	Sphere Power	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Quarterly or Conventional Replacement Lenses										
Hydrasoft® Sphere	methafilcon B / 55%	8.3	14.2	+10.00D to -30.00D (0.25D steps)	4	Daily	18.8	Single and 3-pack vials	No Trials (Singles available for purchase)	
		8.6	15.0				15.7			
Hydrasoft® Sphere Thin	methafilcon B / 55%	8.3	14.2	+10.00D to -30.00D (0.25D steps)	4	Daily	31.3	Single and 3-pack vials	No Trials (Singles available for purchase)	
		8.6	15.0							
Hydrasoft® Aphakic	methafilcon B / 55%	8.3	14.2	+10.25D to +30.00D (0.25D steps)	4	Daily	18.8	Single and 3-pack vials	No Trials (Singles available for purchase)	
		8.6	15.0				15.7			
Hydrasoft® Aphakic Thin	methafilcon B / 55%	8.3	14.2	+10.25D to +30.00D (0.25D steps)	4	Daily	18.8	Single and 3-pack vials	No Trials (Singles available for purchase)	
		8.6	15.0				15.7			



VISIBILITY  
TINT



LOW  
MODULUS



NATURALLY  
WETTABLE



HIGH OXYGEN  
TRANSMISSIBILITY



ASPHERIC  
OPTICS

































OPTIMIZED  
BALLAST DESIGN



UV  
BLOCKING

TORIC LENSES

Brand	Material/ H <sub>2</sub> O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
<b>1-Day Disposables</b>												
clariti® 1 day toric 	somofilcon A / 56%	8.6	14.3	Plano to - 9.00D (0.50D steps after -6.00D)	-0.75, -1.25, -1.75	10°, 20°, 60° 70°, 80°, 90° 100°, 110°, 120° 160°, 170°, 180°	SiHy	Daily Disposable	57	90-pack and 30-pack	1 day toric silicone hydrogel/ 5pk	     
				+0.25D to +4.00D (0.25D steps)		20°, 70°, 90° 110°, 160°, 180°						
				Plano to -9.00D (0.50D steps after -6.00D)	-2.25	10°, 20°, 90° 160°, 170°, 180°						
MyDay® toric 	stenfilcon A / 54%	8.6	14.5	Plano to -10.00D (0.50D steps after -6.00D)	-0.75 -1.25 -1.75	10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°	SiHy	Daily Disposable	80	90-pack	54% 1 day/ Silicone Hydrogel toric/ 5pk	       
+0.50D to +6.00D (0.50D steps)	-2.25	10°, 20°, 90°, 160°, 170°, 180°										
<b>2-Week Replacement Lenses</b>												
Avaira Vitality™ toric 	fanfilcon A / 55%	8.5	14.5	+8.00D to -10.00D (0.50D steps after ± 6.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	SiHy	Daily	90	6-pack	fanfilcon A/ single	     
Biomedics® Toric 	ocufilcon D / 55%	8.7	14.5	+6.00D to -9.00D (0.50D steps after -6.00D and +5.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	4	Daily or Extended	17.8	6-pack	toric contact lens/ single	  
Vertex Toric®  <p><i>Product being discontinued 11/1/19</i></p>	methafilcon A / 55%	8.6	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	4	Daily	17.9	6-pack	No Trials	 



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS











OPTIMIZED BALLAST DESIGN



UV BLOCKING

TORIC LENSES

Brand	Material/ H <sub>2</sub> O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Monthly Replacement Lenses												
Biofinity® toric 	comfilcon A / 48%	8.7	14.5	+8.00D to -10.00D (0.50D steps after +/-6.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	SiHy	Daily or up to 7 days extended	116	6-pack	comfilcon A toric / single	
Biofinity® XR toric 	comfilcon A / 48%	8.7	14.5	+10.00D to +8.50D (0.50D steps) +10.00D to -10.00D (0.50D steps after +/- 6.00D)	-0.75 -1.25 -1.75 -2.25 -2.75 -3.25 -3.75 -4.25 -4.75 -5.25 -5.75	5° to 180° (in 5° steps)	SiHy	Daily or up to 7 days extended	116	6-pack	comfilcon A toric / single	
Proclear® toric 	omafilcon B / 62%	8.8 8.4	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)*	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	2	Daily	30.9	6-pack	Proclear toric/ single	
Proclear® toric XR 	omafilcon B / 62%	8.8 8.4	14.4	+6.25D to +10.00D (0.50D steps after +6.50D) -8.50D to -10.00D (0.50D steps)* +10.00D to -10.00D (0.50D steps after +/-6.50D)*	-0.75 -1.25 -1.75 -2.25 -2.75 to -5.75 (0.50D steps)	5° to 180° (5° steps)	2	Daily	30.9	6-pack	Proclear toric/ single	



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS



OPTIMIZED BALLAST DESIGN



UV BLOCKING

**TORIC LENSES**

Brand	Material/ H <sub>2</sub> O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
<b>Quarterly or Conventional Replacement Lenses</b>												
Hydrasoft® Toric	methafilcon B / 55%	8.3 8.6	14.2	+30.00D to -30.00D (0.25D steps)	-0.50 to -10.00 (0.25D steps)	1° to 180° (1° steps)	4	Daily	14.5	Single and 3-pack vials	No Trials (Singles available for purchase)	
		8.3 8.6 8.9 9.2	15.0									
Hydrasoft® Toric Thin	methafilcon B / 55%	8.3 8.6	14.2	+30.00D to -30.00D (0.25D steps)	-0.50 to -10.00 (0.25D steps)	1° to 180° (1° steps)	4	Daily	14.5	Single and 3-pack vials	No Trials (Singles available for purchase)	
		8.3 8.6 8.9 9.2	15.0									

\* Proclear toric XR is available in the full range of sphere powers and cylinder powers. Sphere powers over +6.00D and -8.00D would be considered XR for all listed cylinder powers. Lenses with a -2.75 to -5.75 cylinder with sphere powers UNDER +6.00D or -8.00D also would be within the XR range.





VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS
















OPTIMIZED BALLAST DESIGN



UV BLOCKING

MULTIFOCAL LENSES

Brand	Material/ H <sub>2</sub> O Content	Base Curve (mm)	Dia (mm)	Sphere Power	ADD Power	Design	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/ Pack Size	Features
<b>1-Day Disposables</b>												
clariti® 1 day multifocal 	somofilcon A / 56%	8.6	14.1	+5.00D to -6.00D (0.25D steps)	Low: up to +2.25D High: +2.50D to +3.00D	center near, progressive intermediate and peripheral distance	SiHy	Daily Disposable	86	90-pack and 30-pack	1 day multifocal silicone hydrogel / 10 pk	    
Proclear® 1 day multifocal 	omafilcon A / 60%	8.7	14.2	+6.00D to -10.00D (0.50D steps after -6.00D)	SINGLE POWER PROFILE	center near aspheric	2	Daily Disposable	28	90-pack	60% 1 Day Proclear multifocal / 10 pk	 
<b>Monthly Replacement Lenses</b>												
Biofinity® multifocal 	comfilcon A / 48%	8.6	14.0	+6.00D to -10.00D (0.50D steps after -6.00D)	+1.00 +1.50 +2.00 +2.50	D lens N lens	SiHy	Daily or up to 7 days extended	142	6-pack	comfilcon A multifocal / single	   
Proclear® multifocal 	omafilcon B / 62%	8.7	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)	+1.00 +1.50 +2.00 +2.50	D lens N lens	2	Daily	21.3	6-pack	Proclear multifocal / single	 



VISIBILITY TINT



LOW MODULUS



NATURALLY WETTABLE



HIGH OXYGEN TRANSMISSIBILITY



ASPHERIC OPTICS






OPTIMIZED BALLAST DESIGN



UV BLOCKING

MULTIFOCAL LENSES

Brand	Material/ H <sub>2</sub> O Content	Base Curve (mm)	Dia (mm)	Sphere Power	ADD Power	Design	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Proclear® multifocal XR  	omafilcon B / 62%	8.7	14.4	+6.25D to +20.00D -8.50D to -20.00D (0.50D steps after +/- 6.50D)	+1.00 +1.50 +2.00 +2.50	D lens N lens	2	Daily	21.3	6-pack	Proclear multifocal / single	
		8.7	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D)	+3.00 +3.50 +4.00							
		8.4	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D)	+1.00, +1.50, +2.00, +2.50, +3.00, +3.50, +4.00							
Proclear® multifocal toric  	omafilcon B / 62%	8.8 8.4	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D) CYL: -0.75 to -5.75 (0.50D steps) AXIS: 5° to 180° (5° steps)	+1.00 to +4.00 (0.50D steps)	D lens N lens	2	Daily	48.6	6-pack	Proclear multifocal toric / single	