

Product Reference Guide

Look to the people, products, and services of CooperVision®.

CONTENTS

Spheres	2 – 3
Torics	4 – 6
Multifocals	7 – 8
Toric Multifocals	8

Customer Service: 800 341 2020
www.coopervision.com



CooperVision®
Live Brightly.®

Product Reference Guide



VISIBILITY
TINT



LOW
MODULUS



NATURALLY
WETTABLE



HIGH OXYGEN
TRANSMISSIBILITY



ASPHERIC
OPTICS



OPTIMIZED
BALLAST DESIGN











UV
BLOCKING



DESIGNED FOR
MYOPIA CONTROL

SPHERE LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
1-Day Disposables										
clariti® 1 day 	somofilcon A / 56%	8.6	14.1	-0.50D to -10.00D +0.50D to +8.00D (0.50D steps after +/- 6.00D) No Plano	SiHy	Daily Disposable	86	90-pack and 30-pack	56% 1 day Silicone Hydrogel sphere / 10 pk	
MyDay® 	stenfilcon A / 54%	8.4	14.2	+8.00D to -12.00D (0.50D steps after +5.00D and -6.00D) No Plano	SiHy	Daily Disposable	100	90-pack and 180-pack	54% 1 day Silicone Hydrogel asphere / 5pk	
Proclear® 1 day 	omafilcon A / 60%	8.7	14.2	+8.00D to -12.00D (0.50D steps after +5.00D and -6.00D) No Plano	2	Daily Disposable	25	90-pack	60% 1 day sphere / 10pk	
MiSight® 1 day 	omafilcon A / 60%	8.7	14.2	-0.50D to -7.00D (0.50D steps after -6.00D) No plano	2	Daily Disposable	25	90-pack	Myopia Management ActivControl™ Technology / 10pk	

Product Reference Guide



VISIBILITY
TINT



LOW
MODULUS



NATURALLY
WETTABLE



HIGH OXYGEN
TRANSMISSIBILITY



ASPHERIC
OPTICS



OPTIMIZED
BALLAST DESIGN





































UV
BLOCKING



DESIGNED FOR
MYOPIA CONTROL

SPHERE LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
2-Week Replacement Lenses										
<div>Avaira Vitality™</div> <div></div>	fanfilcon A / 55%	8.4	14.2	+8.00D to -12.00D (0.50D steps after +/-6.00D) No plano	SiHy	Daily wear	110	6-pack	fanfilcon A / single	<div></div>
<div>Biomedics® 55 Premier asphere</div> <div></div>	ocufilcon D / 55%	8.6 8.9	14.2	-0.25D to -10.00D (0.50D steps after -6.00D) No plano	4	Daily wear	28.1	6-pack	ocufilcon D 55% asphere / 3pk	<div></div>
		8.8	14.2	+0.25D to +8.00D (0.50D steps after +5.00D) No plano						
Monthly Replacement Lenses										
<div>Biofinity®</div> <div></div>	comfilcon A / 48%	8.6	14.0	+8.00D to -12.00D (0.50D steps after ±6.00D) No plano	SiHy	Daily Wear or Extended Wear up to 6 nights/ 7 days; monthly replacement	160	6-pack	comfilcon A / single	<div></div>
<div>Biofinity® XR</div> <div></div>	comfilcon A / 48%	8.6	14.0	+8.50D to +15.00D -12.50D to -20.00D (0.50D steps) No plano	SiHy	Daily Wear or Extended Wear up to 6 nights/ 7 days; monthly replacement	160	6-pack	comfilcon A XR / single	<div></div>
<div>Biofinity Energys®</div> <div></div>	comfilcon A / 48%	8.6	14.0	+8.00D to -12.00D (0.50D steps after ±6.00D) No plano	SiHy	Daily Wear or Extended Wear up to 6 nights/ 7 days; monthly replacement	160	6-pack	Biofinity Energys® / single	<div></div>
<div>Proclear® sphere</div> <div></div>	omafilcon B / 62%	8.6	14.2	+20.00D to -20.00D (0.50D steps after +/-6.00D)	2	Daily wear	36	6-pack	Proclear / 3pk	<div></div>
		8.2	14.2	Plano to -10.00D (0.50D steps after -6.00D)						

Product Reference Guide



VISIBILITY
TINT



LOW
MODULUS



NATURALLY
WETTABLE



HIGH OXYGEN
TRANSMISSIBILITY



ASPHERIC
OPTICS



OPTIMIZED
BALLAST DESIGN

















UV
BLOCKING



DESIGNED FOR
MYOPIA CONTROL

TORIC LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
1-Day Disposables												
clariti® 1 day toric 	somofilcon A / 56%	8.6	14.3	plano to -6.00D (0.25D steps)	-0.75, -1.25, -1.75	Full circle in 10 degrees	SiHy	Daily Disposable	57	90-pack and 30-pack	56% 1 day Silicone Hydrogel toric/ 5pk	     
					-2.25	10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°						
				-6.50D to -9.00D (0.50D steps)	-0.75, -1.25, -1.75	10°, 20°, 60°, 70°, 80°, 90°, 100°, 110°, 120°, 160°, 170°, 180°						
					-2.25	10°, 20°, 90°, 160°, 170°, 180°						
				+0.25D to +4.00D (0.25D steps)	-0.75, -1.25, -1.75	10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°						
MyDay® toric 	stenfilcon A / 54%	8.6	14.5	Plano to -6.00D (0.25D steps)	-0.75 -1.25 -1.75	Full circle in 10° steps	SiHy	Daily Disposable	80	90-pack	54% 1 day/ Silicone Hydrogel toric / 5pk	     
				Plano to -6.00D (0.25D steps)								
				-6.50D to -10.00D (0.50D steps)	-2.25	10°, 20°, 90°, 160°, 170°, 180°						
				+0.50D to +6.00D (0.50D steps)								
				-6.50D to -10.00D (0.50D steps)	-0.75 -1.25 -1.75	10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°						
				+0.25D to +6.00D (0.25D steps)								

Product Reference Guide



VISIBILITY
TINT



LOW
MODULUS



NATURALLY
WETTABLE



HIGH OXYGEN
TRANSMISSIBILITY



ASPHERIC
OPTICS



OPTIMIZED
BALLAST DESIGN














UV
BLOCKING



DESIGNED FOR
MYOPIA CONTROL

TORIC LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
2-Week Replacement Lenses												
Avaira Vitality™ toric 	fanfilcon A / 55%	8.5	14.5	+8.00D to -10.00D (0.50D steps after ± 6.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	SiHy	Daily wear	90	6-pack	fanfilcon A / single	     
Biomedics® Toric 	ocufilcon D / 55%	8.7	14.5	+6.00D to -10.00D (0.50D steps after -6.00D and +5.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	4	Daily wear	17.8	6-pack	ocufilcon D toric / single	  

Product Reference Guide



VISIBILITY
TINT



LOW
MODULUS



NATURALLY
WETTABLE



HIGH OXYGEN
TRANSMISSIBILITY



ASPHERIC
OPTICS



OPTIMIZED
BALLAST DESIGN











UV
BLOCKING



DESIGNED FOR
MYOPIA CONTROL

TORIC LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Monthly Replacement Lenses												
Biofinity® toric 	comfilcon A / 48%	8.7	14.5	+8.00D to -10.00D (0.50D steps after +/-6.00D)	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	SiHy	Daily Wear or Extended Wear up to 6 nights/ 7 days; monthly replacement	116	6-pack	comfilcon A toric / single	
Biofinity® XR toric 	comfilcon A / 48%	8.7	14.5	+8.50D to +20.00D and -10.50D to -20.00D (0.50D steps) +20.00D to -20.00D (0.50D steps after +/-6.00D)	-0.75 -1.25 -1.75 -2.25 -2.75 -3.25 -3.75 -4.25 -4.75 -5.25 -5.75	5° to 180° (in 5° steps)	SiHy	Daily Wear or Extended Wear up to 6 nights/ 7 days; monthly replacement	116	6-pack	comfilcon A toric / single	
Proclear® toric 	omafilcon B / 62%	8.8 8.4	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)*	-0.75 -1.25 -1.75 -2.25	10° to 180° (10° steps)	2	Daily wear	25	6-pack	Proclear® toric / single	
Proclear® toric XR 	omafilcon B / 62%	8.8 8.4	14.4	+6.25D to +10.00D (0.50D steps after +6.50D) -8.50D to -10.00D (0.50D steps)* +10.00D to -10.00D (0.50D steps after +/-6.50D)*	-0.75 -1.25 -1.75 -2.25 -2.75 to -5.75 (0.50D steps)	5° to 180° (5° steps)	2	Daily wear	39	6-pack	Proclear® toric / single	

Product Reference Guide



VISIBILITY
TINT



LOW
MODULUS



NATURALLY
WETTABLE



HIGH OXYGEN
TRANSMISSIBILITY



ASPHERIC
OPTICS



OPTIMIZED
BALLAST DESIGN













UV
BLOCKING



DESIGNED FOR
MYOPIA CONTROL

MULTIFOCAL LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	ADD Power	Design	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
1-Day Disposables												
clariti® 1 day multifocal 	somofilcon A / 56%	8.6	14.1	+6.00D to -8.00D in 0.25 D steps	Low: up to +2.25D High: +2.50D to +3.00D	center near, progressive intermediate and peripheral distance	SiHy	Daily Disposable	86	90-pack and 30-pack	56% 1 day Silicone Hydrogel multifocal / 10 pk	
MyDay® multifocal 	stenfilcon A / 54%	8.4	14.2	+8.00D to -10.00D (0.25 D steps) -10.50D to -12.00D (0.50 D steps)	Low: +0.75D to +1.25D Mid: +1.50D to +1.75D High: +2.00D to +2.50D	Binocular Progressive System™	SiHy	Daily Disposable	100	90-pack	54% 1 day Silicone Hydrogel multifocal / 10pk	
Proclear® 1 day multifocal 	omafilcon A / 60%	8.7	14.2	+6.00D to -10.00D (0.50D steps after -6.00D)	SINGLE POWER PROFILE	center near aspheric	2	Daily Disposable	28	90-pack	60% 1 day PC multifocal/ 10 pk	
Monthly Replacement Lenses												
Biofinity® multifocal 	comfilcon A / 48%	8.6	14.0	+6.00D to -10.00D (0.50D steps after -6.00D)	+1.00 +1.50 +2.00 +2.50	D lens N lens	SiHy	Daily Wear or Extended Wear up to 6 nights/ 7 days; monthly replacement	142	6-pack	comfilcon A multifocal / single	
Proclear® multifocal 	omafilcon B / 62%	8.7	14.4	+6.00D to -8.00D (0.50D steps after -6.50D)	+1.00 +1.50 +2.00 +2.50	D lens N lens	2	Daily wear	17	6-pack	Proclear® multifocal / single	

Product Reference Guide



VISIBILITY
TINT



LOW
MODULUS



NATURALLY
WETTABLE



HIGH OXYGEN
TRANSMISSIBILITY



ASPHERIC
OPTICS



OPTIMIZED
BALLAST DESIGN






UV
BLOCKING










DESIGNED FOR
MYOPIA CONTROL

MULTIFOCAL LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	ADD Power	Design	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Monthly Replacement Lenses												
Proclear® multifocal XR 	omafilcon B / 62%	8.7	14.4	+6.25D to +20.00D -8.50D to -20.00D (0.50D steps after +/- 6.50D)	+1.00 +1.50 +2.00 +2.50	D lens N lens	2	Daily wear	17	6-pack	Proclear multifocal / single	 
		8.7	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D)	+3.00 +3.50 +4.00							
		8.4	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D)	+1.00, +1.50, +2.00, +2.50, +3.00, +3.50, +4.00							

TORIC MULTIFOCAL LENSES

Brand	Material/ H ₂ O Content	Base Curve (mm)	Dia (mm)	Sphere Power	Cylinder Power	Axis	ADD Power	Design	FDA Group	Wearing Schedule	Dk/t (at -3.00D)	Revenue Carton Size	Trials Labeled As/Pack Size	Features
Monthly Replacement Lenses														
Biofinity® toric multifocal 	comfilcon A / 48%	8.7	14.5	+10.00D to -10.00D (0.50D steps after +/- 6.00D)	-0.75D to -5.75D (0.50D steps)	5°-180°, 5° steps	+1.00 +1.50 +2.00 +2.50	D lens N lens	SiHy	Daily Wear or Extended Wear up to 6 nights/7 days; monthly replacement	116	6-pack	Biofinity® toric multifocal / single	    
Proclear® multifocal toric 	omafilcon B / 62%	8.8 8.4	14.4	+20.00D to -20.00D (0.50D steps after +/- 6.50D)	-0.75 to -5.75 (0.50D steps)	5° to 180° (5° steps)	+1.00 to +4.00 (0.50D steps)	D lens N lens	2	Daily wear	38.6	6-pack	Proclear® multifocal toric / single	