



CooperVision®
Live Brightly.®

Product Specifications Guide



MyDay® daily disposable family

MyDay® daily disposable



Wearing Schedule

Daily disposable

Sphere Power

+8.00D to -12.00D
(0.50D steps after
+5.00D and
-6.00D)

No plano

Carton Size Options

90 pack
180 pack

Design

Aspheric

Material / H₂O Content

stenfilcon A / 54%

Base Curve 8.4mm

Diameter 14.2mm

FDA Group SiHy

UV Blocking* Yes

Visibility Tint Yes

High-performing 1-day contact lenses for demanding lifestyles.

MyDay® daily disposable toric



Wearing Schedule

Daily disposable

Sphere Power

Plano to -10.00D
(0.50D steps after
-6.00D)

+0.50D to +6.00D
(0.50D steps)

Carton Size Options

90 pack

Cylinder Power

-0.75
-1.25
-1.75

-2.25

Material / H₂O Content

stenfilcon A / 54%

Axis

10°, 20°, 70°,
80°, 90°, 100°,
110°, 160°, 170°,
180°

10°, 20°, 90°,
160°, 170°, 180°

Base Curve 8.6mm

Diameter 14.5mm

FDA Group SiHy

UV Blocking* Yes

Visibility Tint Yes

Design

Optimized
Toric Lens
Geometry™

High-performing 1-day contact lenses for wearers with astigmatism.

clariti® 1 day family

clariti® 1 day



Wearing Schedule

Daily disposable

Sphere Power

-0.50D to -10.00D
+0.50D to +8.00D
(0.50D steps after
±6.00D)

No plano

Carton Size Options

90 pack
30 pack

Design

Aspheric

Material / H₂O Content

somofilcon A / 56%

Base Curve 8.6mm

Diameter 14.1mm

FDA Group SiHy

UV Blocking* Yes

Visibility Tint No

Great everyday value and high breathability¹ in a convenient daily disposable contact lens.

clariti® 1 day toric



Great everyday value and high breathability¹ in a convenient toric daily disposable contact lens.

Wearing Schedule Daily disposable	Carton Size Options 90 pack 30 pack	Material / H₂O Content somofilcon A / 56%	Base Curve 8.6mm Diameter 14.3mm FDA Group SiHy UV Blocking* Yes Visibility Tint No Design Aspheric
Sphere Power Plano to -6.00D (0.25D steps)	Cylinder Power -0.75 -1.25 -1.75 -2.25	Axis 10° to 180° (10° steps) 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°	
-6.50D to -9.00D (0.50D steps)	-0.75 -1.25 -1.75 -2.25	10°, 20°, 60°, 70°, 80°, 90°, 100°, 110°, 120°, 160°, 170°, 180° 10°, 20°, 90°, 160°, 170°, 180°	
+0.25D to +4.00D (0.25D steps)	-0.75 -1.25 -1.75	10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°	

clariti® 1 day multifocal



Great everyday value and high breathability¹ in the convenience of a daily disposable multifocal contact lens.

Wearing Schedule Daily disposable	Carton Size Options 90 pack 30 pack	Material / H₂O Content somofilcon A / 56%	Base Curve 8.6mm Diameter 14.1mm FDA Group SiHy UV Blocking* Yes Visibility Tint No
Sphere Power +5.00D to -6.00D (0.25D steps)	Add Power Low: up to +2.25D High: +2.50D to +3.00D	Design Center near, progressive intermediate, and peripheral distance	

Biofinity® family

Biofinity Energys®



Wearing Schedule

Daily or 6 nights/
7 days Extended
wear

Sphere Power

+8.00D to -12.00D
(0.50D steps after
±6.00D)

Carton Size Options

6 pack

Design

Aspheric
Digital Zone
Optics® lens
design

Material / H₂O Content

comfilcon A / 48%

Base Curve 8.6mm

Diameter 14.0mm

FDA Group SiHy

UV Blocking* No

Visibility Tint Yes

Contact lenses for your digital
lifestyle.

Biofinity®



Wearing Schedule

Daily or 6 nights/
7 days Extended
wear

Sphere Power

+15.00D to
-20.00D
(0.50D steps after
+/-6.00D)

Carton Size Options

6 pack

Design

Asphere

Material / H₂O Content

comfilcon A / 48%

Base Curve 8.6mm

Diameter 14.0mm


FDA Group SiHy

UV Blocking* No

Visibility Tint Yes

Remarkable visual performance.
Premium all-day comfort.

Biofinity® toric



Wearing Schedule

Daily or up to 6
nights/7 days
Extended wear

Sphere Power

+8.00D to -10.00D
(0.50D steps after
+/-6.00D)

Carton Size Options

6 pack

Cylinder Power

-0.75
-1.25
-1.75
-2.25

Material / H₂O Content

comfilcon A / 48%

Axis

10° to 180°
(in 10° steps)

Base Curve 8.7mm

Diameter 14.5mm

FDA Group SiHy

UV Blocking* No

Visibility Tint Yes

Design
Optimized
Toric Lens
Geometry™

Biofinity® XR toric



+8.50D to +20.00D
and -10.50D
to -20.00D
(0.50D steps)

-0.75
-1.25
-1.75
-2.25

5° to 180°
(in 5° steps)


+20.00D to
-20.00D
(0.50D steps after
+/-6.00D)

-2.75
-3.25
-3.75
-4.25
-4.75
-5.25
-5.75

5° to 180°
(in 5° steps)

Remarkable visual performance.
Premium all-day comfort.

Biofinity[®] multifocal



Wearing Schedule

Daily wear or
Extended wear up
to 6 nights/7 days

Sphere Power

+6.00D to -10.00D
(0.50D steps after
-6.00D)

Carton Size Options

6 pack

Add Power

+1.00
+1.50
+2.00
+2.50

Material / H₂O Content

comfilcon A / 48%

Design

D lens
N lens

Base Curve 8.6mm
Diameter 14.0mm
FDA Group SiHy
UV Blocking* No
Visibility Tint Yes

Remarkable visual performance at every distance. Premium all-day comfort.

Biofinity[®] toric multifocal



Wearing Schedule

Daily wear or
Extended wear up
to 6 nights/7 days;
Monthly
replacement

Sphere Power

+10.00D to
-10.00D (0.50D
steps after
+/-6.00D)

Carton Size Options

6 pack

Add Power

+1.00
+1.50
+2.00
+2.50

Material / H₂O Content

comfilcon A / 48%

Design


D lens
N lens

Base Curve 8.7mm
Diameter 14.5mm
FDA Group SiHy
UV Blocking* No
Visibility Tint Yes

Delivering performance to contact lens wearers who are both astigmatic and presbyopic.

Avaira Vitality[™] family

Avaira Vitality[™]



Wearing Schedule

Daily

Sphere Power

+8.00D to -12.00D
(0.50D steps after
+/-6.00D)

Carton Size Options

6 pack

Design

Aspheric front
surface

Material / H₂O Content


fanfilcon A / 55%

Base Curve 8.4mm
Diameter 14.2mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint Yes

High performance. Outstanding value.

No plano

Avaira Vitality™ toric



Wearing Schedule
Daily

Carton Size Options
6 pack

Material / H₂O Content
fanfilcon A / 55%

Base Curve 8.5mm

Diameter 14.5mm

FDA Group SiHy

UV Blocking* Yes

Visibility Tint Yes

Sphere Power
+8.00D to -10.00D
(0.50D steps after
±6.00D)

Cylinder Power
-0.75
-1.25
-1.75
-2.25

Axis
10° to 180°
(10° steps)

Design
Optimized Toric
Lens Geometry™

High performance. Outstanding value.

Brilliant Futures™ Myopia Management Program with MiSight® 1 day

Brilliant Futures™ Myopia Management Program with MiSight®



Wearing Schedule
Daily

Carton Size Options
90 pack

Material / H₂O Content
omafilcon A / 60%

Base Curve 8.7mm

Diameter 14.2mm

FDA Group 2

UV Blocking* No


Visibility Tint Yes

Sphere Power
-0.50D to -6.00D
(0.25D steps)

Child-friendly daily disposable contact lenses proven to slow the progression of nearsightedness in children aged 8-12 at the initiation of treatment.^{2†}

Proclear® family

Proclear® 1 day



Wearing Schedule

Daily disposable

Carton Size Options

90 pack

Material / H₂O Content

omafilcon A / 60%

Base Curve 8.7mm

Diameter 14.2mm

FDA Group 2

UV Blocking* No

Visibility Tint Yes

Sphere Power

+8.00D to -12.00D
(0.50D steps
after -6.00D and
+5.00D)


Design

Aspheric

No plano

Clear vision³, convenience and comfort, for contact lens wearers who experience dry eyes due to lens wear.

Proclear® 1 day multifocal



Wearing Schedule

Daily disposable

Carton Size Options

90 pack

Material / H₂O Content

omafilcon A / 60%

Base Curve 8.7mm

Diameter 14.2mm

FDA Group 2

UV Blocking* No

Visibility Tint Yes

Sphere Power

+6.00D to -10.00D
(0.50D steps after
-6.00D)

Add Power

Single power
profile

Design

Center near
aspheric

Convenience and comfort for multifocal contact lens wearers who experience dry eyes due to lens wear.

Proclear®



Wearing Schedule

Daily

Carton Size Options

6 pack

Material / H₂O Content

omafilcon B / 62%

Base Curve 8.6, 8.2mm

Diameter 14.2mm

FDA Group 2

UV Blocking* No

Visibility Tint Yes

Sphere Power

+20.00D to
-20.00D (0.5D
steps after +/-
6.00D)

Design

Sphere

Clear vision³ and comfort for contact lens wearers who experience dry eyes due to lens wear.

Plano to -10.00D
(0.50D steps after
-6.00D)

Proclear® toric



Wearing Schedule
Daily

Carton Size Options
6 pack

Material / H₂O Content
omafilcon B / 62%

Base Curve 8.8, 8.4mm

Diameter 14.4mm

FDA Group 2

UV Blocking* No

Visibility Tint Yes

Sphere Power
+6.00D to -8.00D
(0.50D steps after
-6.50D)

Cylinder Power
-0.75
-1.25
-1.75
-2.25

Axis
10° to 180°
(10° steps)

Helps provide all-day comfort for contact lens wearers with astigmatism who experience dry eyes due to lens wear.

Proclear® multifocal



Wearing Schedule
Daily

Carton Size Options
6 pack

Material / H₂O Content
omafilcon B / 62%

Base Curve 8.7mm

Diameter 14.4mm

FDA Group 2

UV Blocking* No

Visibility Tint Yes

Sphere Power
+6.00D to -8.00D
(0.50D steps after
-6.50D)

Add Power
+1.00
+1.50
+2.00
+2.50

Design
D lens
N lens

Helps provide all-day comfort for multifocal contact lens wearers who experience dry eyes related to lens wear.

Proclear® multifocal toric



Wearing Schedule
Daily

Carton Size Options
6 pack

Material / H₂O Content
omafilcon B / 62%

Base Curve 8.8, 8.4mm

Diameter 14.4mm

FDA Group 2

UV Blocking* No

Visibility Tint Yes

Sphere Power
+20.00D to
-20.00D (0.50D
steps after +/-
6.50D)

Cylinder Power
-0.75 to -5.75
(0.50D steps)

Axis
5° to 180°
(5° steps)

Design
D lens
N lens

Helps provide all-day comfort for multifocal contact lens wearers with astigmatism who experience dry eyes related to lens wear.

Biomedics® family

Biomedics® 55 Premier



Designed for comfort and performance at an affordable price.

Wearing Schedule
Daily or Extended
(1 to 7 days
between removals)

Sphere Power
-0.25D to -10.00D
(0.50D steps after
-6.00D)

+0.25D to +6.00D
(0.50D steps after
+5.00D)


Carton Size Options
6 pack

Design
Aspheric

Material / H₂O Content
ocufilcon D / 55%

Base Curve 8.6, 8.9, 8.8mm
Diameter 14.2mm
FDA Group 4
UV Blocking* No
Visibility Tint Yes

Biomedics® toric



A quality contact lens for wearers with astigmatism.

Wearing Schedule
Daily or Extended
(1 to 7 days
between removals)

Sphere Power
+6.00D to -9.00D
(0.50D steps
after -6.00D and
+5.00D)

Carton Size Options
6 pack

Cylinder Power
-0.75
-1.25
-1.75
-2.25

Material / H₂O Content
ocufilcon D / 55%

Axis
10° to 180°
(10° steps)

Base Curve 8.7mm
Diameter 14.5mm
FDA Group 4
UV Blocking* Yes
Visibility Tint Yes

¹ Brennan NA. Beyond Flux: Total Corneal Oxygen Consumption as an Index of Corneal Oxygen During Contact Lens Wear. *Optom Vis Sci* 2005. 2005;82(6):467-472.

² Chamberlain P, et al. A 3-year randomized clinical trial of MiSight® lenses for myopia control. *Optom Vis Sci*. 2019;96(8):556-567.

³ CVI data on file 2017 & 2014; non-dispensing short-term fitting study. N=40 & n=20.

* Compared with a single-vision one-day lens over a three-year period.

*UV-blocking contact lenses help provide protection against transmission of harmful UV radiation to eye but are not substitutes for protective UV-absorbing eyewear, such as UV absorbing goggles or sunglasses, as they do not completely cover the eye or surrounding area. Continue to use UV-absorbing eyewear as directed by your eye care professional.