



CooperVision®  
Live Brightly.®

# Product Specifications Guide



## MyDay® daily disposable family

### MyDay® daily disposable


**Wearing Schedule**

Daily disposable

**Sphere Power**

+8.00D to -12.00D  
(0.50D steps after  
+5.00D and  
-6.00D)

No plano

**Carton Size Options**

90 pack  
180 pack

**Design**

Aspheric

**Material / H<sub>2</sub>O Content**

stenfilcon A / 54%

**Base Curve** 8.4mm

**Diameter** 14.2mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** Yes

High-performing 1-day contact lenses for demanding lifestyles.

### MyDay® daily disposable toric


**Wearing Schedule**

Daily disposable

**Sphere Power**

Plano to -10.00D  
(0.50D steps after  
-6.00D)

+0.50D to +6.00D  
(0.50D steps)

**Carton Size Options**

90 pack

**Cylinder Power**

-0.75  
-1.25  
-1.75

-2.25

**Material / H<sub>2</sub>O Content**

stenfilcon A / 54%

**Axis**

10°, 20°, 70°,  
80°, 90°, 100°,  
110°, 160°, 170°,  
180°

10°, 20°, 90°,  
160°, 170°, 180°

**Base Curve** 8.6mm

**Diameter** 14.5mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** Yes

**Design**

Optimized  
Toric Lens  
Geometry™

High-performing 1-day contact lenses for wearers with astigmatism.

## clariti® 1 day family

### clariti® 1 day


**Wearing Schedule**

Daily disposable

**Sphere Power**

-0.50D to -10.00D  
+0.50D to +8.00D  
(0.50D steps after  
±6.00D)

No plano

**Carton Size Options**

90 pack  
30 pack

**Design**

Aspheric

**Material / H<sub>2</sub>O Content**

somofilcon A / 56%

**Base Curve** 8.6mm

**Diameter** 14.1mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** No

Great everyday value and high breathability<sup>1</sup> in a convenient daily disposable contact lens.

## clariti® 1 day toric



Great everyday value and high breathability<sup>1</sup> in a convenient toric daily disposable contact lens.

<b>Wearing Schedule</b>	<b>Carton Size Options</b>	<b>Material / H<sub>2</sub>O Content</b>	<b>Base Curve</b>	8.6mm
Daily disposable	90 pack 30 pack	somofilcon A / 56%	<b>Diameter</b>	14.3mm
<b>Sphere Power</b>	<b>Cylinder Power</b>	<b>Axis</b>	<b>FDA Group</b>	SiHy
Plano to -6.00D (0.25D steps)	-0.75 -1.25 -1.75  -2.25	10° to 180° (10° steps)  10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°	<b>UV Blocking*</b>	Yes
-6.50D to -9.00D (0.50D steps)	-0.75 -1.25 -1.75  -2.25	10°, 20°, 60°, 70°, 80°, 90°, 100°, 110°, 120°, 160°, 170°, 180°  10°, 20°, 90°, 160°, 170°, 180°	<b>Visibility Tint</b>	No
+0.25D to +4.00D (0.25D steps)	-0.75 -1.25 -1.75	10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°	<b>Design</b>	Aspheric

## clariti® 1 day multifocal



Great everyday value and high breathability<sup>1</sup> in the convenience of a daily disposable multifocal contact lens.

<b>Wearing Schedule</b>	<b>Carton Size Options</b>	<b>Material / H<sub>2</sub>O Content</b>	<b>Base Curve</b>	8.6mm
Daily disposable	90 pack 30 pack	somofilcon A / 56%	<b>Diameter</b>	14.1mm
<b>Sphere Power</b>	<b>Add Power</b>	<b>Design</b>	<b>FDA Group</b>	SiHy
+5.00D to -6.00D (0.25D steps)	Low: up to +2.25D  High: +2.50D to +3.00D	Center near, progressive intermediate, and peripheral distance	<b>UV Blocking*</b>	Yes
			<b>Visibility Tint</b>	No

## Biofinity® family

### Biofinity Energys®



#### Wearing Schedule

Daily or 6 nights/  
7 days Extended  
wear

#### Sphere Power

+8.00D to -12.00D  
(0.50D steps after  
±6.00D)

#### Carton Size Options

6 pack

#### Design

Aspheric  
Digital Zone  
Optics® lens  
design

#### Material / H<sub>2</sub>O Content

comfilcon A / 48%

**Base Curve** 8.6mm

**Diameter** 14.0mm

**FDA Group** SiHy

**UV Blocking\*** No

**Visibility Tint** Yes

Contact lenses for your digital  
lifestyle.

### Biofinity®



#### Wearing Schedule

Daily or 6 nights/  
7 days Extended  
wear

#### Sphere Power

+15.00D to  
-20.00D  
(0.50D steps after  
+/-6.00D)

#### Carton Size Options

6 pack

#### Design

Asphere

#### Material / H<sub>2</sub>O Content

comfilcon A / 48%

**Base Curve** 8.6mm

**Diameter** 14.0mm

**FDA Group** SiHy

**UV Blocking\*** No

**Visibility Tint** Yes

Remarkable visual performance.  
Premium all-day comfort.

### Biofinity® toric



#### Wearing Schedule

Daily or up to 6  
nights/7 days  
Extended wear

#### Sphere Power

+8.00D to -10.00D  
(0.50D steps after  
+/-6.00D)

#### Carton Size Options

6 pack

#### Cylinder Power

-0.75  
-1.25  
-1.75  
-2.25

#### Material / H<sub>2</sub>O Content

comfilcon A / 48%

#### Axis

10° to 180°  
(10° steps)

**Base Curve** 8.7mm

**Diameter** 14.5mm

**FDA Group** SiHy

**UV Blocking\*** No

**Visibility Tint** Yes

#### Design

Optimized  
Toric Lens  
Geometry™

Remarkable visual performance.  
Premium all-day comfort.

+10.00D to +8.50D  
(0.50D steps)

-0.75  
-1.25  
-1.75  
-2.25

5° to 180°  
(5° steps)

+10.00D to  
-10.00D  
(0.50D steps after  
+/- 6.00D)

-2.75  
-3.25  
-3.75  
-4.25  
-4.75  
-5.25  
-5.75

5° to 180°  
(5° steps)

## Biofinity® multifocal



### Wearing Schedule

Daily wear or  
Extended wear up  
to 6 nights/7 days

### Sphere Power

+6.00D to -10.00D  
(0.50D steps after  
-6.00D)

### Carton Size Options

6 pack

### Add Power

+1.00  
+1.50  
+2.00  
+2.50

### Material / H<sub>2</sub>O Content

comfilcon A / 48%

### Design

D lens  
N lens

**Base Curve** 8.6mm  
**Diameter** 14.0mm  
**FDA Group** SiHy  
**UV Blocking\*** No  
**Visibility Tint** Yes

Remarkable visual performance at  
every distance. Premium all-day  
comfort.

## Biofinity® toric multifocal



### Wearing Schedule

Daily wear or  
Extended wear up  
to 6 nights/7 days;  
Monthly  
replacement

### Sphere Power

+10.00D to  
-10.00D (0.50D  
steps after  
+/-6.00D)

### Carton Size Options

6 pack

### Add Power

+1.00  
+1.50  
+2.00  
+2.50

### Material / H<sub>2</sub>O Content

comfilcon A / 48%

### Design


D lens  
N lens

**Base Curve** 8.7mm  
**Diameter** 14.5mm  
**FDA Group** SiHy  
**UV Blocking\*** No  
**Visibility Tint** Yes

Delivering performance to contact  
lens wearers who are both  
astigmatic and presbyopic.

## Avaira Vitality™ family

### Avaira Vitality™



### Wearing Schedule

Daily

### Sphere Power

+8.00D to -12.00D  
(0.50D steps after  
+/-6.00D)

### Carton Size Options

6 pack

### Design

Aspheric front  
surface

### Material / H<sub>2</sub>O Content

fanfilcon A / 55%

**Base Curve** 8.4mm  
**Diameter** 14.2mm  
**FDA Group** SiHy  
**UV Blocking\*** Yes  
**Visibility Tint** Yes

High performance. Outstanding  
value.

## Avaira Vitality™ toric



**Wearing Schedule**  
Daily

**Sphere Power**  
+8.00D to -10.00D  
(0.50D steps after  
±6.00D)

**Carton Size Options**  
6 pack

**Cylinder Power**  
-0.75  
-1.25  
-1.75  
-2.25

**Material / H<sub>2</sub>O Content**  
fanfilcon A / 55%

**Axis**  
10° to 180°  
(10° steps)

**Design**  
Optimized Toric  
Lens Geometry™

**Base Curve** 8.5mm  
**Diameter** 14.5mm  
**FDA Group** SiHy  
**UV Blocking\*** Yes  
**Visibility Tint** Yes

High performance. Outstanding value.

## Brilliant Futures™ with MiSight® Myopia Management

### Brilliant Futures™ with MiSight® Myopia Management



**Wearing Schedule**  
Daily

**Sphere Power**  
-6.00D to -0.25D  
(0.25D steps)

**Carton Size Options**  
30 pack

**Material / H<sub>2</sub>O Content**  
omafilcon A

**Base Curve** 8.7mm  
**Diameter** 14.2mm  
**FDA Group** 2  
**UV Blocking\*** Yes  
**Visibility Tint** Yes

Child-friendly daily disposable contact lenses proven to slow the progression of nearsightedness in children aged 8-12 at the initiation of treatment.<sup>2‡</sup>

## Proclear® family

### Proclear® 1 day



#### Wearing Schedule

Daily disposable

#### Carton Size Options

90 pack

#### Material / H<sub>2</sub>O Content

omafilcon A / 60%

**Base Curve** 8.7mm

**Diameter** 14.2mm

**FDA Group** 2

**UV Blocking\*** No

**Visibility Tint** Yes

#### Sphere Power

+8.00D to -12.00D  
(0.50D steps  
after -6.00D and  
+5.00D)


#### Design

Aspheric

No plano

Clear vision<sup>3</sup>, convenience and comfort, for contact lens wearers who experience dry eyes due to lens wear.

### Proclear® 1 day multifocal



#### Wearing Schedule

Daily disposable

#### Carton Size Options

90 pack

#### Material / H<sub>2</sub>O Content

omafilcon A / 60%

**Base Curve** 8.7mm

**Diameter** 14.2mm

**FDA Group** 2

**UV Blocking\*** No

**Visibility Tint** Yes

#### Sphere Power

+6.00D to -10.00D  
(0.50D steps after  
-6.00D)

#### Add Power

Single power  
profile

#### Design

Center near  
aspheric

Convenience and comfort for multifocal contact lens wearers who experience dry eyes due to lens wear.

### Proclear®



#### Wearing Schedule

Daily

#### Carton Size Options

6 pack

#### Material / H<sub>2</sub>O Content

omafilcon B / 62%

**Base Curve** 8.6, 8.2mm

**Diameter** 14.2mm

**FDA Group** 2

**UV Blocking\*** No

**Visibility Tint** Yes

#### Sphere Power

+20.00D to  
-20.00D (0.5D  
steps after +/-  
6.00D)

#### Design

Sphere

Clear vision<sup>3</sup> and comfort for contact lens wearers who experience dry eyes due to lens wear.

Plano to -10.00D  
(0.50D steps after  
-6.00D)

## Proclear® toric



**Wearing Schedule**  
Daily

**Carton Size Options**  
6 pack

**Material / H<sub>2</sub>O Content**  
omafilcon B / 62%

**Base Curve** 8.8, 8.4mm

**Diameter** 14.4mm

**FDA Group** 2

**UV Blocking\*** No

**Visibility Tint** Yes

**Sphere Power**  
+6.00D to -8.00D  
(0.50D steps after  
-6.50D)

**Cylinder Power**  
-0.75  
-1.25  
-1.75  
-2.25

**Axis**  
10° to 180°  
(10° steps)

Helps provide all-day comfort for contact lens wearers with astigmatism who experience dry eyes due to lens wear.

## Proclear® multifocal



**Wearing Schedule**  
Daily

**Carton Size Options**  
6 pack

**Material / H<sub>2</sub>O Content**  
omafilcon B / 62%

**Base Curve** 8.7mm

**Diameter** 14.4mm

**FDA Group** 2

**UV Blocking\*** No

**Visibility Tint** Yes

**Sphere Power**  
+6.00D to -8.00D  
(0.50D steps after  
-6.50D)

**Add Power**  
+1.00  
+1.50  
+2.00  
+2.50

**Design**  
D lens  
N lens

Helps provide all-day comfort for multifocal contact lens wearers who experience dry eyes related to lens wear.

## Proclear® multifocal toric



**Wearing Schedule**  
Daily

**Carton Size Options**  
6 pack

**Material / H<sub>2</sub>O Content**  
omafilcon B / 62%

**Base Curve** 8.8, 8.4mm

**Diameter** 14.4mm

**FDA Group** 2

**UV Blocking\*** No

**Visibility Tint** Yes

**Sphere Power**  
+20.00D to  
-20.00D (0.50D  
steps after +/-  
6.50D)

**Cylinder Power**  
-0.75 to -5.75  
(0.50D steps)

**Axis**  
5° to 180°  
(5° steps)

Helps provide all-day comfort for multifocal contact lens wearers with astigmatism who experience dry eyes related to lens wear.

**Design**  
D lens  
N lens



## Biomedics<sup>®</sup> family

### Biomedics<sup>®</sup> 55 Premier



Designed for comfort and performance at an affordable price.

**Wearing Schedule**  
Daily or Extended  
(1 to 7 days  
between removals)

**Sphere Power**  
-0.25D to -10.00D  
(0.50D steps after  
-6.00D)  
  
+0.25D to +6.00D  
(0.50D steps after  
+5.00D)

**Carton Size  
Options**  
6 pack

**Design**  
Aspheric

**Material / H<sub>2</sub>O  
Content**  
ocufilcon D / 55%

**Base Curve** 8.6, 8.9, 8.8mm  
**Diameter** 14.2mm  
**FDA Group** 4  
**UV Blocking\*** No  
**Visibility Tint** Yes

### Biomedics<sup>®</sup> toric



A quality contact lens for wearers  
with astigmatism.

**Wearing Schedule**  
Daily or Extended  
(1 to 7 days  
between removals)

**Sphere Power**  
+6.00D to -9.00D  
(0.50D steps  
after -6.00D and  
+5.00D)

**Carton Size  
Options**  
6 pack

**Cylinder  
Power**  
-0.75  
-1.25  
-1.75  
-2.25

**Material / H<sub>2</sub>O  
Content**  
ocufilcon D / 55%

**Axis**  
10° to 180°  
(10° steps)

**Base Curve** 8.7mm  
**Diameter** 14.5mm  
**FDA Group** 4  
**UV Blocking\*** Yes  
**Visibility Tint** Yes

<sup>1</sup> Brennan NA. Beyond Flux: Total Corneal Oxygen Consumption as an Index of Corneal Oxygen During Contact Lens Wear. Optom Vis Sci 2005. 2005;82(6):467-472.

<sup>2</sup> Chamberlain P, et al. A 3-year randomized clinical trial of MiSight<sup>®</sup> lenses for myopia control. Optom Vis Sci. 2019;96(8):556-567.

<sup>3</sup> CVI data on file 2017 & 2014; non-dispensing short-term fitting study. N=40 & n=20.

\* Compared with a single-vision one-day lens over a three-year period.

\*UV-blocking contact lenses help provide protection against transmission of harmful UV radiation to eye but are not substitutes for protective UV-absorbing eyewear, such as UV absorbing goggles or sunglasses, as they do not completely cover the eye or surrounding area. Continue to use UV-absorbing eyewear as directed by your eye care professional.

All trademarks, service marks, trade names, trade dress and product names featured herein are registered trademarks and trademarks of CooperVision, its affiliates and subsidiaries, unless otherwise noted. CooperVision's trademarks may be used publicly only with permission.