



## EASILY FIT MORE PRESBYOPES MORE PRECISELY. NOW AND FOR YEARS TO COME.

Frequency's heritage of exceptional visual acuity continues with the excellent performance of Frequency 55 Multifocal. With four ADD powers, Frequency 55 Multifocal lenses allow you to fit a wider range of presbyopes, from early to advanced, and to keep them in Frequency 55 Multifocal lenses as their condition progresses.

**BALANCED PROGRESSIVE™ TECHNOLOGY** — two different, complementary lenses work together to produce exceptional vision at every distance.

**UNMATCHED FITTING FLEXIBILITY** — four ADD powers and a wide range of sphere powers let you fit patients more precisely.

**CONSISTENT COMFORT** — patented UltraSync® manufacturing provides superior lens reproducibility and a thin, comfortable edge.

**SUPERIOR PATIENT RETENTION** — availability of higher ADD powers allows you to continue fitting patients as their presbyopia progresses.







## THE POWER OF BALANCED PROGRESSIVE TECHNOLOGY.

Balanced Progressive Technology combines spherical and aspheric optics and unique zone sizes to produce a “D” lens for the dominant eye and an “N” lens for the non-dominant eye. These two lens designs work together to provide a new level of vision quality.

### D LENS: DOMINANT EYE



### N LENS: NON-DOMINANT EYE



## PROFESSIONAL-GRADE CONTACT LENSES.

Frequency 55 Multifocal lenses' advanced Balanced Progressive Technology allows each lens to be fit independently. Combine that with a greater number of ADD powers, and only Frequency 55 Multifocal offers the level of fitting flexibility and visual acuity that uncompromising professionals like you demand.

Put it all together, and you'll see that Frequency 55 Multifocal lenses are a worthy addition to the renowned Frequency family of lenses, each one offering your patients unsurpassed visual quality.

PARAMETERS	MATERIAL	WEARING SCHEDULE	WATER CONTENT	BASE CURVE (mm)	DIAMETER (mm)	MFG. PROCESS	SPHERE POWER	ADD POWER	DESIGN	HANDLING TINT	LENS CARE
<b>FREQUENCY 55 MULTIFOCAL</b>	methafilcon A	Flexible wear	55%	8.7	14.4	UltraSync	+4.00D to -6.00D (in 0.25D steps)	+1.00 +1.50 +2.00 +2.50	D lens N lens	Light blue	Chemical or peroxide disinfection

**FOR A MULTIFOCAL FITTING CONSULTATION, CALL THE CVSPO MULTIFOCAL FITTING TEAM AT 866-316-9358 OR VISIT COOPERVISION.COM/MULTIFOCAL.**